



Toxicology

One simple device for saliva drug testing.

# iScreen<sup>®</sup> OFD

Oral Fluid Drug Screening Device.



For forensic use only.

Substance abuse testing with more **substance**.

## How it works.

PRODUCT TRAINING AND  
CERTIFICATION AVAILABLE ONLINE.

### Features & benefits

- » Accurately screen for illicit and prescription drugs in saliva
- » Test remains in control of administrator
- » No bathroom necessary
- » Simple collection and testing procedure
- » Results in as little as 10 minutes
- » Long shelf life at room temperature storage
- » Integrated internal procedural control

Remove Sample Collector Protector and insert Sponge end of collector into mouth.




Actively swab the inside of the mouth and tongue with Sponge.



Insert the collector vertically into the Cap and twist the Handle clockwise 180° to tighten the Cap. Read results at 10 minutes.



If results are presumptive positive, secure Cap with Security Seal and send to laboratory for confirmation.



PRODUCT PROCEDURE NOTE: Refer to product insert for complete instructions, limitations, and warnings.

This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS), gas chromatography/tandem mass spectrometry (GC/MS/MS) and liquid chromatography/tandem mass spectrometry (LC/MS/MS) are the preferred confirmatory methods. Clinical consideration and professional judgement should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

For forensic use only.

### Oral Fluid Screening Device Cutoff Levels

Test	Calibrator	Cutoff
Amphetamine (AMP)	d-Amphetamine	50 ng/mL
Cocaine (COC)	Benzoylcegonine	20 ng/mL
Marijuana (THC)	11-nor- $\Delta^9$ -THC-9-COOH	12 ng/mL
Methamphetamine (mAMP)	d-Methamphetamine	50 ng/mL
Opiate (OPI)	Morphine	40 ng/mL
Phencyclidine (PCP)	Phencyclidine	10 ng/mL

