










Orawell® Oral Fluid Drug Screen Device Reference Guide

The Orawell® Oral Fluid Drug Screen Device is a simple one-step test for the detection of drugs of abuse in oral fluid.

<p>1. Remove the blue cap by pulling the both sides gently. This will expose the collection pad.</p> 	<p>2. Gently rub the collection pad inside mouth against cheek in a circular motion several (approximately 10-15) times.</p> 	<p>3. Gently rub the collection pad against the opposite cheek in circular motion several (approximately 10-15) times.</p> 
<p>4. Gently rub the collection pad on top of the tongue and underneath of the tongue several (approximately 10-15) times. Do not chew, suck, bite or bend the collection pad.</p> 	<p>5. Place the collection pad underneath the tongue for approximately 30 seconds to collect saliva. Instruct the donor to hold the device in place with hand. When sufficient amount of saliva is collected, the red fluid will move up to the window.</p> 	<p>6. Remove from mouth after red fluid shows up in both of the test windows. Re-cap the device and lay on a flat surface. Start timer and read results in 5 minutes.</p> 

Interpreting Test Results

<p style="text-align: center;"><u>Invalid Result</u></p> <p>When no colored band appears in the CONTROL (C) region, the test is invalid even if there is a band in the test region. Repeat the test with a new device.</p> <p>Example Interpretation:</p>  <table style="margin-left: 20px;"> <tr><td>COC:</td><td>Invalid</td></tr> <tr><td>OPI:</td><td>Invalid</td></tr> <tr><td>MET:</td><td>Invalid</td></tr> <tr><td>BZO:</td><td>Invalid</td></tr> <tr><td>THC:</td><td>Invalid</td></tr> <tr><td>AMP:</td><td>Invalid</td></tr> </table>	COC:	Invalid	OPI:	Invalid	MET:	Invalid	BZO:	Invalid	THC:	Invalid	AMP:	Invalid	<p style="text-align: center;"><u>Negative Result</u></p> <p>For each test, two colored bands should be observed:</p> <ul style="list-style-type: none"> • One in the CONTROL (C) region • One in the specific TEST region <p>The color of the test band may be slightly darker or lighter than the control band. Any visible band that can be seen is a negative result.</p> <p>Example Interpretation:</p>  <table style="margin-left: 20px;"> <tr><td>COC:</td><td>Negative</td></tr> <tr><td>OPI:</td><td>Negative</td></tr> <tr><td>MET:</td><td>Negative</td></tr> <tr><td>BZO:</td><td>Negative</td></tr> <tr><td>THC:</td><td>Negative</td></tr> <tr><td>AMP:</td><td>Negative</td></tr> </table>	COC:	Negative	OPI:	Negative	MET:	Negative	BZO:	Negative	THC:	Negative	AMP:	Negative	<p style="text-align: center;"><u>Presumptive Positive Result</u></p> <p>A colored band at the CONTROL (C) region should be observed. When there is no colored band at the specific TEST region, the test is presumptive positive for that particular drug.</p> <p>Example Interpretation:</p>  <table style="margin-left: 20px;"> <tr><td>COC:</td><td>Negative</td></tr> <tr><td>OPI:</td><td>Negative</td></tr> <tr><td>MET:</td><td>Presumptive Positive</td></tr> <tr><td>BZO:</td><td>Negative</td></tr> <tr><td>THC:</td><td>Negative</td></tr> <tr><td>AMP:</td><td>Negative</td></tr> </table>	COC:	Negative	OPI:	Negative	MET:	Presumptive Positive	BZO:	Negative	THC:	Negative	AMP:	Negative
COC:	Invalid																																					
OPI:	Invalid																																					
MET:	Invalid																																					
BZO:	Invalid																																					
THC:	Invalid																																					
AMP:	Invalid																																					
COC:	Negative																																					
OPI:	Negative																																					
MET:	Negative																																					
BZO:	Negative																																					
THC:	Negative																																					
AMP:	Negative																																					
COC:	Negative																																					
OPI:	Negative																																					
MET:	Presumptive Positive																																					
BZO:	Negative																																					
THC:	Negative																																					
AMP:	Negative																																					