


Dependable results and simple to administer.

Alere iCassette[®]

Urine Drug Screening Device



For professional in vitro diagnostic use only


Knowing now matters.TM

Alere iCassette® Urine Drug Screen

Features and Benefits

- Rapidly screen up to 11 illicit drugs and prescription drugs
- Administrator remains in control of test device
- Simple procedure: collect and pipette to test—results in as little as 5 minutes

Test in 3 Easy Steps



PRODUCT PROCEDURE NOTE: Refer to product insert for complete instructions, limitations, and warnings. For healthcare professionals including professionals at point of care sites. For in vitro diagnostic use only. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography-mass spectrometry (GC-MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

Drug screen options and cutoff levels

Test	Calibrator	Cutoff Level
Amphetamine (AMP)	d-Amphetamine	1,000 ng/mL
Barbiturates (BAR)	Secobarbital	300 ng/mL
Benzodiazepines (BZO)	Oxazepam	300 ng/mL
Cocaine (COC)	Benzoyllecgonine	300 ng/mL
Marijuana (THC)	11-nor- Δ^9 -THC-9 COOH	50 ng/mL
Methadone (MTD)	Methadone	300 ng/mL
Methamphetamine (mAMP)	d-Methamphetamine	1,000 ng/mL
Methylenedioxymethamphetamine (MDMA - Ecstasy)	d,l Methylenedioxymethamphetamine	500 ng/mL
Morphine (MOP 300 or OPI 300)	Morphine	300 ng/mL
Opiates (OPI 2000)	Morphine	2,000 ng/mL
Phencyclidine (PCP)	Phencyclidine	25 ng/mL
Tricyclic Antidepressants (TCA)	Nortriptyline	1,000 ng/mL



Products 800.340.4029 | aleretoxicology.com